



**pediatric**  
**BIOBRANE® GLOVES**  
Temporary Wound Dressing

The Biobrane® Glove is intended for covering clean partial thickness burn wounds of the hand. The glove promotes early movement of the burned hand. **CAUTION:** The Biobrane® Glove is not recommended for meshed autografts on fingers or hands because of the potential for shearing the grafted tissue.

Please see Glove Instructions available on this website.



Distributed by  
**UDL LABORATORIES INC.**  
A subsidiary of Mylan Laboratories Inc.